



Support arm – fixed

Support arm R1100000

Support arm with leg, use height 850 mm R1101000


Support arm with leg, use height 800 mm R1102000

Provides firm support when you are sitting, standing, or transferring to and from a toilet. Good, broad support surface for the hand and arm. The rounded form gives a firm and comfortable grip. The support arm can fold up when not in use or to make more space.




Mounting: R1100: When it is mounted at the recommended height (890/840 (DK) mm), the support arm's height is 850/800 mm – measured from the floor to the top edge of the support arm.

Fixing bolts for solid walls are included.



Mounting: R1101/1102: When it is mounted at the recommended height (890/840 (DK) mm), the support arm's height is 850/800 mm – measured from the floor to the top edge of the support arm.

Fixing bolts for solid walls are included.



Support arms for mounting on toilet


With toilet seat R37 R1170

Without toilet seat R1171

Can be mounted easily and quickly on the toilet. The folding arms give a good support and the moulded handles are comfortable for hands and arms.


The support arms can be folded individually.

White toilet seat and cobalt handles.  
Supplied with hinges BL1.



Accessory products for R1100/R1101/R1102


R1151



Paper holder: With spring. The resistance of the spring enables it to be used with one hand.

Accessory products for R1170/R1171

R1175



Paper holder: Can be mounted at the right or the left handle.

Rxxxxx



Toilet seat  
See page 20.

- 000 White
- Grey, end of handle